


# Get Started with Suki on mobile


Copy/Paste process prior to Epic Integration

**\*Obtain verbal consent from the patient prior to recording\***

## 1 View your schedule



From the home page, tap Select a patient or tap Patients to view a list of patients you've already seen.




Enter the patient's first and last name. Tap +New Patient to add a patient that doesn't appear in your list.

## 2 Start ambient




Tap Start Ambient to begin recording the visit.



Ambient recording is started. Tap Pause if you need to leave the room or have an unrelated discussion.


## 3 Conduct the visit




### For Improved Notes:

- Set mobile device down between you and the patient
- Talk through physical exam and diagnoses
- Summarize your plan

## 4 End ambient



Tap the check mark to stop recording the visit.



Choose a note type to help Suki determine the sections of content to generate.